

Product Datasheet

Rabbit Anti-Goat IgG (H+L) Secondary Antibody, Cy3 Conjugated (orb27706)

Description	Rabbit Anti-Goat IgG (H+L) Secondary Antibody, Cy3 Conjugated
Species/Host	Rabbit
Conjugation	Cy3
Tested Applications	FC, IF
Concentration	1mg/ml
Storage	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.