

Product Datasheet

Rabbit Anti-Goat IgG (H+L) Secondary Antibody, Cy3 Conjugated (orb27706)

Description	Cy3 Conjugated Rabbit Anti-Goat IgG (H+L) secondary antibody This Cy3 conjugated antibody is specific for goat IgG and shows no cross-reactivity with human/rat/mouse/rabbit IgG.
Species/Host	Rabbit
Conjugation	Cy3
Tested Applications	FC, IF
Concentration	1mg/ml
Storage	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)