

## Product Datasheet

# Rabbit IgG Isotype Control (H+L) Antibody (orb27739)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Rabbit IgG Isotype Control (H+L)                           |
| <b>Species/Host</b>    | Rabbit   |
| <b>Conjugation</b>     | Unconjugated   |
| <b>Storage</b>         | At -20°C for 2 years. Avoid repeated freezing and thawing. |
| <b>Note</b>            | For research use only                                      |
| <b>Clonality</b>       | Polyclonal   |
| <b>Expiration Date</b> | 12 months from date of receipt.                            |