


## Product Datasheet


# Polyclonal Antibody against Human Lipocalin-2 (LCN2) antibody (orb27770)

<b>Description</b>	Rabbit polyclonal antibody to polyclonal
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.




P-Peptide -

Immunohistochemical analysis of paraffin...



P-Peptide +

Immunohistochemical analysis of paraffin...



Western blot analysis of extract from He...