

Product Datasheet

CD79A antibody (orb27792)

Description	Mouse monoclonal antibody to CD79A is found in the development from the earliest pre-B cell to plasm...
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CD79A
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	25038
Uniprot ID	P11912
NCBI	NP_067612.1 , NP_001774.1
Dilution Range	WB: 1:100-250
Expiration Date	12 months from date of receipt.

