



Product Datasheet

CNOT4 antibody (orb27802)

Description	Mouse monoclonal antibody to CNOT4 has an E3 ubiquitin ligase activity that is involved in activatio...
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	CNOT4
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG2a,k
Clonality	Monoclonal
MW	63510
Uniprot ID	O95628
NCBI	NP_001177779.1 , NP_001177778.1 , NP_001177777.1 , NP_001177776.1 , NP_001008226.1 , NP_037448.2
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse heart tis...