

Product Datasheet

Ribosome-inactivating protein PD-L1/PD-L2 antibody (orb27820)

Description	Ribosome-inactivating protein PD-L1/PD-L2 antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Phytolacca dioica Ribosome-inactivating protein PD-L1/PD-L2 protein (1-261AA)
Target	Ribosome-inactivating protein PD-L1/PD-L2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P84853
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.