

Product Datasheet

Amatoxin antibody (orb27830)

Description	Amatoxin antibody
Species/Host	Rabbit
Reactivity	Fungi
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Peptide sequence from Amanita phalloides Alpha-amanitin proprotein protein
Target	Alpha-amanitin proprotein
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P85421
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.