


Product Datasheet

Ribosome-inactivating protein cucurmosin antibody (orb27831)

Description	Ribosome-inactivating protein cucurmosin antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Cucurbita moschata Ribosome-inactivating protein cucurmosin protein (1-244AA)
Target	Ribosome-inactivating protein cucurmosin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P84531
Dilution Range	WB:1:1000-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Cucurbita mosch...