

Product Datasheet

Ribosome-inactivating protein cucurmosin antibody (HRP) (orb27832)

| | |
|----------------------------|---|
| Description | Ribosome-inactivating protein cucurmosin antibody (HRP) |
| Species/Host | Rabbit |
| Reactivity | Plant |
| Conjugation | HRP |
| Tested Applications | ELISA |
| Immunogen | Recombinant Cucurbita moschata Ribosome-inactivating protein cucurmosin protein (1-244AA) |
| Target | Ribosome-inactivating protein cucurmosin |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P84531 |
| Dilution Range | ELISA |
| Expiration Date | 12 months from date of receipt. |