

## Product Datasheet

# Ribosome-inactivating protein cucurmosin antibody (Biotin) (orb27834)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Ribosome-inactivating protein cucurmosin antibody (Biotin)                                |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Plant   |
| <b>Conjugation</b>         | Biotin  |
| <b>Tested Applications</b> | ELISA   |
| <b>Immunogen</b>           | Recombinant Cucurbita moschata Ribosome-inactivating protein cucurmosin protein (1-244AA) |
| <b>Target</b>              | Ribosome-inactivating protein cucurmosin  |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol                            |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                             |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Uniprot ID</b>          | <a href="#">P84531</a>  |
| <b>Dilution Range</b>      | ELISA   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |