

Product Datasheet

Major pollen allergen Dac g 4 antibody (orb27840)

| | |
|----------------------------|---|
| Description | Major pollen allergen Dac g 4 antibody |
| Species/Host | Rabbit |
| Reactivity | Plant |
| Conjugation | Unconjugated |
| Tested Applications | ELISA |
| Immunogen | Recombinant Dactylis glomerata Major pollen allergen Dac g 4 protein (1-55AA) |
| Target | Major pollen allergen Dac g 4 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P82946 |
| Dilution Range | ELISA |
| Expiration Date | 12 months from date of receipt. |