

Product Datasheet

Major pollen allergen Dac g 4 antibody (HRP) (orb27842)

Description	Major pollen allergen Dac g 4 antibody (HRP)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Dactylis glomerata Major pollen allergen Dac g 4 protein (1-55AA)
Target	Major pollen allergen Dac g 4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P82946
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.