

Product Datasheet

Pollen allergen Dac g 3 antibody (Biotin) (orb27868)

Description	Pollen allergen Dac g 3 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Dactylis glomerata Pollen allergen Dac g 3 protein (1-96AA)
Target	Pollen allergen Dac g 3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P93124
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.