

Product Datasheet

Ribosome-inactivating protein karasurin-C antibody (FITC) (orb27976)

| | |
|------------------------|--|
| Description | Ribosome-inactivating protein karasurin-C antibody (FITC) |
| Species/Host | Rabbit |
| Reactivity | Plant |
| Conjugation | FITC |
| Immunogen | Recombinant Trichosanthes kirilowii Ribosome-inactivating protein karasurin-C protein (22-270AA) |
| Target | Ribosome-inactivating protein karasurin-C |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P24478 |
| Dilution Range | ELISA |
| Expiration Date | 12 months from date of receipt. |