



Product Datasheet

DYKDDDDK Tag antibody (orb27997)

Description	Rabbit polyclonal antibody to V5 Tag is derived from a small epitope, which presents on the P and V ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	KLH conjugated synthetic peptide (CDYKDDDDKDYKDDDDK) was used as antigen.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



WB analysis of 12 tag protein using DYKD...