

Product Datasheet

PF-06647263 Antibody (orb2803553)

Description	PF-06647263 (anti-EFNA4-ADC) is an ADC compound targeting EFNA4, with anti-tumor activity, composed of an anti-EFNA4 antibody conjugated with calicheamicin, inducing significant tumor regression in TNBC xenografts.
Conjugation	Unconjugated
Target	PF-06647263
Concentration	4.56 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Purity	> 95%
MW	145.5 kDa
CAS Number	1822383-65-2
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)