

Product Datasheet

Protein-C Tag antibody (orb28041)

Description	Rabbit polyclonal antibody to NNMT belongs to the NNMT/PNMT/TEMT family that is highly expressed in ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PROC
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	52071
Uniprot ID	P04070
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

