

Product Datasheet

Mouse Apelin (APLN) ELISA Kit (orb2809777)

Description	Mouse Apelin (APLN) ELISA Kit
Reactivity	Mouse
Range	125-8000pg/mL
Note	For research use only
Application notes	Assay Type: Competitive
Uniprot ID	Q9R0R4
Sensitivity	64pg/mL
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)