

Product Datasheet

Ribosome-inactivating protein bryodin I antibody (orb28117)

Description	Ribosome-inactivating protein bryodin I antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Bryonia dioica Ribosome-inactivating protein bryodin I protein (24-247AA)
Target	Ribosome-inactivating protein bryodin I
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P33185
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.