

Product Datasheet

Ribosome-inactivating protein bryodin I antibody (Biotin) (orb28120)

Description	Ribosome-inactivating protein bryodin I antibody (Biotin)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Bryonia dioica Ribosome-inactivating protein bryodin I protein (24-247AA)
Target	Ribosome-inactivating protein bryodin I
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P33185
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.