

## Product Datasheet

### Human Interleukin 35 (IL-35) ELISA Kit (orb2813082)

|                          |                                        |
|--------------------------|----------------------------------------|
| <b>Description</b>       | Human Interleukin 35 (IL-35) ELISA Kit |
| <b>Reactivity</b>        | Human                                  |
| <b>Range</b>             | 31.25-2000pg/mL                        |
| <b>Note</b>              | For research use only                  |
| <b>Application notes</b> | Assay Type: Sandwich                   |
| <b>Sensitivity</b>       | 16.25pg/mL                             |
| <b>Expiration Date</b>   | Please enquire.                        |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)