

Product Datasheet

Thalidomide-piperidine-C2-piperazine-Boc (orb2814494)

Description	Thalidomide-piperidine-C2-piperazine-Boc is an E3 Ligase Ligand-linker conjugate comprising a cereblon ligand based on Thalidomide and a linker. This compound is utilized in the synthesis of PROTAC.
Storage	-20°C
Note	For research use only
Formula	$C_{29}H_{39}N_5O_6$
MW	553.65
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)