

## Product Datasheet

### GNAS antibody (HRP) (orb28211)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | GNAS antibody (HRP)  |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human  |
| <b>Conjugation</b>         | HRP  |
| <b>Tested Applications</b> | ELISA  |
| <b>Immunogen</b>           | Recombinant Human Guanine nucleotide-binding protein G(s) subunit alpha isoforms short protein (1-394AA) |
| <b>Target</b>              | GNAS   |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol  |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | <a href="#">P63092</a>   |
| <b>Dilution Range</b>      | ELISA  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |