

Product Datasheet

Ribosome-inactivating protein alpha-trichosanthin antibody (Biotin) (orb28301)

Description	Ribosome-inactivating protein alpha-trichosanthin antibody (Biotin)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Trichosanthes kirilowii Ribosome-inactivating protein alpha-trichosanthin protein (24-270AA)
Target	Ribosome-inactivating protein alpha-trichosanthin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P09989
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.