




Product Datasheet

CF150 antibody (orb28416)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	CGAS
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	58814
Uniprot ID	Q8N884
NCBI	NP_612450.2
Dilution Range	FACS: 1:10-50, WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of MDA-MB231 ce...



WB analysis of MCF-7 HT-29 THP-1 cell li...