

Product Datasheet

Trypsin inhibitor DE-3 antibody (FITC) (orb28465)

Description	Trypsin inhibitor DE-3 antibody (FITC)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	FITC
Immunogen	Recombinant Erythrina caffra Trypsin inhibitor DE-3 protein (1-172AA)
Target	Trypsin inhibitor DE-3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P09943
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.