

Product Datasheet

ACF antibody (orb28588)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications FC, IHC-P, WB

Immunogen Synthetic Peptide

Target A1CF

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal



MW 65202

Uniprot ID [Q9NQ94](#)


NCBI [NP_001185749.1](#), [NP_001185747.1](#), [NP_055391.2](#), [NP_620311.1](#), [NP_620310.1](#), [NP_001185748.1](#)

Dilution Range WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50

Expiration Date 12 months from date of receipt.



Western blot analysis of A375 cell line ...



Immunohistochemical analysis of formalin...