




Product Datasheet

HIST1H2BJ antibody (orb28652)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	H2BC11
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	13904
Uniprot ID	P06899
NCBI	NP_066402.2
Dilution Range	IHC-P: 1:50-100, FACS: 1:10-50, WB: 1:1000
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of CEM cells (r...



Western blot analysis of CEM, HL-60, NCI...