




Product Datasheet

SLC16A11 antibody (orb28857)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	SLC16A11
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	47791
Uniprot ID	Q8NCK7
NCBI	NP_699188.1
Dilution Range	FACS: 1:10-50, IHC-P: 1:50-100, WB: 1:1000
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of WiDr cells (...)



Western blot analysis of mouse brain and...