

Product Datasheet

GRIN2A antibody (orb29358)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications FC, IHC-P, WB

Immunogen Synthetic Peptide

Target GRIN2A

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal



MW 165283

Uniprot ID [Q12879](#)


NCBI [NP_000824.1](#), [NP_001127879.1](#), [NP_001127880.1](#)

Dilution Range WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50

Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of formalin...



Flow cytometric analysis of MDA-MB435 ce...