

Product Datasheet

SMTNL2 antibody (orb29689)

Description	Rabbit polyclonal antibody to ID1 is an inhibitor of differentiation protein, and is located at the ...
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SMTNL2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	50196
Uniprot ID	Q2TAL5
NCBI	NP_940903.2 , NP_001108446.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

