

## Product Datasheet

### TAG-72 antibody (orb302237)

<b>Description</b>	Mouse monoclonal antibody to TAG-72
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P
<b>Immunogen</b>	Please contact us for more information.
<b>Preservatives</b>	Tris Buffer, pH 7.3-7.7, containing 1% rAlbumin and 0.1% Sodium Azide
<b>Concentration</b>	n/a
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1-k
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	IHC-P (1:100-1:500)
<b>Expiration Date</b>	12 months from date of receipt.