

Product Datasheet

CD32 antibody (FITC) (orb302316)

Description	Mouse monoclonal antibody against CD32 conjugated to FITC
Species/Host	Mouse
Reactivity	Canine, Human, Porcine, Primate
Conjugation	FITC
Tested Applications	FC
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.09% sodium azide and 1% rAlbumin
Concentration	0.1 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.