

Product Datasheet

c-Myc Antibody (orb303233)

Description	Chicken polyclonal antibody to c-Myc
Species/Host	Gallus
Conjugation	Unconjugated
Immunogen	KLH-conjugated peptide, sequence MEQKLISEEDLNE, human c-Myc UniProt: Q6LBK7
Form/Appearance	Lyophilized
Storage	Store lyophilized antibodies at -20°C and reconstituted antibodies at 4°C. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube.
Note	For research use only
Application notes	c-Myc is derived from Myc proto-oncogene protein. A short peptide located between amino acids 408-420 is used as an epitope tag, that is easily recognized by tag specific antibodies. This is a universal detection reagent which will allow detection of any tag containing protein.
Clonality	Polyclonal
Purity	Immunogen affinity purified IgY in PBS pH 8.5. Contains 0.02% azide.
MW	Depends on target protein which is tagged with this epitope
Uniprot ID	P01106
Dilution Range	IF/ICC: 1:2000-5000
Expiration Date	12 months from date of receipt.