




Product Datasheet

POMT1 antibody (orb30384)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	POMT1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	84881
Uniprot ID	Q9Y6A1
NCBI	NP_001070834.1 , NP_009102.3 , NP_001070833.1 , NP_001129586.1 , NP_001129585.1
Dilution Range	IHC-P: 1:10-50, WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of ZR-75-1, K562 c...



Flow cytometric analysis of U-251MG cell...