



## Product Datasheet

### DARPP32 (phospho-T75) antibody (orb304516)

<b>Description</b>	Rabbit polyclonal antibody to DARPP32 (phospho-T75)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Porcine, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IH, WB
<b>Immunogen</b>	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T75 of human DARPP32 protein. The exact sequence is proprietary.
<b>Target</b>	PPP1R1B
<b>Preservatives</b>	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
<b>Form/Appearance</b>	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Source</b>	Rabbit
<b>Uniprot ID</b>	<a href="#">Q9UD71</a> , <a href="#">Q60829</a> , <a href="#">Q6J4I0</a>
<b>Entrez</b>	<a href="#">84152</a> , <a href="#">360616</a> , <a href="#">19049</a>
<b>Dilution Range</b>	WB: 1:500-1000, IHC-P: 1:100-200
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of DARPP32 (Phospho-T75)



Immunohistochemical analysis of DARPP32