

Product Datasheet

NF-kappaB p65 (phospho-T254) antibody (orb304547)

Description

Rabbit polyclonal antibody to NF-kappaB p65



Western blot analysis of NF-kappaB p65 (...

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Primate, Rat

Conjugation

Unconjugated

Tested Applications

IH, IP, WB

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T254 of human NF-kappaB p65 protein. The exact sequence is proprietary.

Target

RELA

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source


Rabbit

Uniprot ID
[Q04207](#), [Q04206](#)
Entrez
[5970](#), [19697](#)
Dilution Range


WB: 1:500-1000, IHC-P: 1:100-200, IP: 1:10-100

Expiration Date

12 months from date of receipt.



Immunohistochemical analysis of NF-kappa...



Direct ELISA antibody dose-response curv...