

Product Datasheet

DARPP32 (phospho-T75) peptide (orb304757)

Description	The peptide is used to DARPP32 (phospho-T75) antibody ractivity
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	PPP1R1B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9UD71 , Q60829 , Q6J4I0
Entrez	84152 , 360616 , 19049
Expiration Date	12 months from date of receipt.