

Product Datasheet

PKC mu (pY463) peptide (orb304762)

Description	The peptide is used to PKC mu (pY463) antibody ractivity
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	PRKD1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9WTQ1 , Q62101 , Q15139
Entrez	85421 , 18760 , 5587
Expiration Date	12 months from date of receipt.