

Product Datasheet

PTP1B (phospho-S50) peptide (orb304771)

Description	The peptide is used to PTP1B (phospho-S50) antibody ractivity
Reactivity	Human, Mouse, Porcine, Primate, Rat, Zebrafish
Tested Applications	Blocking
Target	PTPN1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P18031 , P35821 , P20417
Entrez	24697 , 5770 , 19246
Expiration Date	12 months from date of receipt.