

Product Datasheet

Ryanodine Receptor 2 (phospho-S2808) peptide (orb304815)

Description	The peptide is used to Ryanodine Receptor 2 (phospho-S2808) antibody ractivity
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	RYR2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	E9Q401 , B0LPN4 , Q92736
Entrez	20191 , 6262 , 689560
Expiration Date	12 months from date of receipt.