

Product Datasheet

SHCA (phospho-Y349) peptide (orb304825)

Description	The peptide is used to SHCA (phospho-Y349) antibody ractivity
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SHC1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P98083 , Q5M824 , P29353
Entrez	6464 , 20416 , 85385
Expiration Date	12 months from date of receipt.