

Product Datasheet

SMAD2 (phospho-S467) peptide (orb304834)

Description	The peptide is used to SMAD2 (phospho-S467) antibody ractivity
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	SMAD2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	O70436 , Q62432 , Q15796
Entrez	29357 , 17126 , 4087
Expiration Date	12 months from date of receipt.