

Product Datasheet

SRF (phospho-S103) peptide (orb304848)

Description	The peptide is used to SRF (phospho-S103) antibody ractivity
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	SRF
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9JM73 , P11831
Entrez	6722 , 20807
Expiration Date	12 months from date of receipt.