

Product Datasheet

MST1/2 (phospho-T183) peptide (orb304853)

Description	The peptide is used to MST1/2 (phospho-T183) antibody ractivity
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Tested Applications	Blocking
Target	STK3; STK4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9JI10 , O54748 , Q13188 , Q13043 , Q9JI11
Entrez	58231 , 56274 , 6789 , 65189
Expiration Date	12 months from date of receipt.