

## Product Datasheet

### TCF4/12 peptide (orb304866)

<b>Description</b>	The peptide is used to TCF4/12 antibody ractivity
<b>Reactivity</b>	Human, Mouse, Porcine, Rat
<b>Tested Applications</b>	Blocking
<b>Target</b>	TCF4; TCF12
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">P15884</a> , <a href="#">Q62655</a> , <a href="#">Q61286</a> , <a href="#">P51514</a> , <a href="#">Q60722</a> , <a href="#">Q99081</a>
<b>Entrez</b>	<a href="#">21413</a> , <a href="#">6925</a> , <a href="#">6938</a> , <a href="#">84382</a> , <a href="#">25720</a> , <a href="#">21406</a>
<b>Expiration Date</b>	12 months from date of receipt.