

## Product Datasheet

### ACK1 (phospho-Y284) peptide (orb304878)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | The peptide is used to ACK1 (phospho-Y284) antibody ractivity   |
| <b>Reactivity</b>          | Human, Mouse, Porcine, Rat  |
| <b>Tested Applications</b> | Blocking  |
| <b>Target</b>              | TNK2  |
| <b>Form/Appearance</b>     | Lyophilized powder  |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.  |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours. |
| <b>Purity</b>              | > 85%   |
| <b>Source</b>              | Synthetic   |
| <b>Uniprot ID</b>          | <a href="#">Q07912</a> , <a href="#">O54967</a> , <a href="#">Q5U2X5</a>  |
| <b>Entrez</b>              | <a href="#">10188</a> , <a href="#">303882</a> , <a href="#">51789</a>  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |