

Product Datasheet

p63 (phospho-S455) peptide (orb304886)

Description	The peptide is used to p63 (phospho-S455) antibody ractivity
Reactivity	Human
Tested Applications	Blocking
Target	TP63
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9H3D4
Entrez	8626
Expiration Date	12 months from date of receipt.