

Product Datasheet

UBA5 peptide (orb304895)

Description	The peptide is used to UBA5 antibody ractivity
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	UBA5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9GZZ9 , Q5M7A4 , Q8VE47
Entrez	79876 , 300968 , 66663
Expiration Date	12 months from date of receipt.